



NON STERILE LATEX GLOVES POWDER FREE & POWDERED

Designed for everyday use, our non-sterile gloves offer superior comfort, reliable protection, and a secure fit to keep your hands safe throughout the day.

GLOVES SIZES

Ex-Small, Small, Medium, Large, Ex-Large

QUALITY CONTROL

- Free from dirt marks, oil stains, embedded foreign particles, coagulum etc.
- Ambidextrous with smooth surface fit either hand and textured in the fingertips for better grip.
- The cuff fit well and stays in place without rolling or folding.
- Powder free gloves with a special blend of polymer coating for easy wearing on dry or wet hands.
- Powdered gloves: powder content < 10 mg/dm²
Powder-free gloves: powder residue < 2 mg per glove.

QUALITY STANDARDS

- Conforms to ASTM D3577, EN 455(part-1,2 , 3 and 4) , ISO 10282 , IS 13422 , EN ISO 374 (part 1,2,4 and 5), EN 16523-1 and EN ISO 21420 standards.
- Manufactured under EN ISO 13485:2016,ISO 14001:2015 and ISO 9001:2015 Quality Quality Management System.
- CE certified as per PPE Category III and MDR Class I.

PACKING SPECIFICATIONS

Packing	Quantity
Dispenser Box	100 Pieces
Carton Box	2000 Pieces

Caution ⚠: This product contains natural rubber latex which may cause allergic reactions including anaphylactic responses in some sensitised individuals.



Our Corporate Office and Factory

St. Marys Rubbers Pvt. Ltd XVII/401A, Thottamkavala, Vizhikkathode, Koovappally P.O, Kanjirappally, Kottayam - 686518, Kerala, India. E-mail:info@stmarysrubbers.com, Website:www.stmarysrubbers.com, Customer Care : +91 9447486270, +914828 252277